

CISME Data Analysis Script Walkthrough

1. Download R studio: https://www.rstudio.com/products/rstudio/download/#download




- 2. Make sure CISME output is saved as a .csv file into the same folder as provided R scripts (CISME PR Analysis.R and CISME PR DoNotEdit.Rmd)
- 3. Open the CISME PR Analysis.R script in R studio
- 4. Set the working directory (the location where your CISME data output file and *CISME PR DoNotEdit.Rmd* are found)
 - > Select 'Session' menu
 - > Select 'Set Working Directory'
 - > Select 'To Source File Location'



5. Input CISME output file name (without .csv extension) within the "quotations" at the top of the *CISME PR Analysis*. script.

- 6. Select 'Source' in the top right screen.
- 7. A web browser will pop up to input parameters for analysing the output. The only REQUIRED entry is 'Input File Name' (again, input file name without .csv extension).



- 8. Other parameter options:
 - > Salinity: you can modify salinity here if the salinity in the CISME metadata is not correct
 - ➤ Dark time limit (min): time in minutes at which the user would like to use as a cut off for dark (respiration) analyses
 - ➤ **Light minimum time (min):** time in minutes at which the user would like to use as the start for light (photosynthesis) analyses
 - ➤ **Light maximum time (min):** time in minutes at which the user would like to use as a cut off for light (photosynthesis) analyses

- ➤ Initial measured total alkalinity (uMol/kg): measured initial total alkalinity taken before sample run. * Without total alkalinity measurements, no calculation of TA, total carbon dioxide, or calcification rates will be calculated.
- Final measured total alkalinity (uMol/kg): measured final total alkalinity from CISME sample tube after sample run
- Measured surface area (cm²): used for standardization of calcification rates
- ➤ **Measured biomass:** used for standardization of calcification rates if not using surface area of sample
- ➤ Total water volume: volume (g) of water in the sample loop and CISME device used to calculate calcification rate
- ➤ Dark calcification percent (decimal): user-determined percentage to calculate dark and light portions of calcification rates (if calculating rates from combined photosynthesis/respiration run)
- 9. Once desired parameters are input in web browser, press 'Save' or hit enter. Pull back up R studio window and wait until the console displays:
 - "Output created: your filename here.html"
- 10. A new folder (named: *yourfilename_files_currentdate*) should now appear in the directory where the original CISME output file is saved. The folder will contain the following:
 - \rightarrow HTML overview document (*filename.html*) \rightarrow this file will open in a web browser
 - > Spreadsheet of the regression table (filename RegressionTable currentdate.csv)
 - ➤ PDF of all figures (filename figures currentdate.pdf)
 - Final spreadsheet with calculated parameters (*filename workup currentdate.csv*)